



M2M -Auth



The registration form is titled "Create Your Account" with a subtitle "Join us today and get started". It features four input fields: "First Name" (with a person icon and the text "John"), "Last Name" (with a person icon and the text "Doe"), "Email" (with an envelope icon and the text "john@example.com"), and "Password" (with a lock icon and masked text "*****"). Below the fields is an orange "Register Now" button with a right-pointing arrow. At the bottom, there is a link "Already have an account? Sign in". A small circular logo with the letter "N" is located in the bottom left corner of the page.

SwapDrill Registration Page



A user-friendly and responsive form designed to onboard users with a seamless account creation experience.



The sign-in form is titled "Sign In". It contains two input fields: "Email" (with the placeholder text "Enter your email") and "Password" (with the placeholder text "Enter your password"). Below the fields is an orange "Sign In" button. At the bottom, there is a link "Dont have an account? Sign up".


SwapDrill Sign-In Page

A secure and intuitive login interface enabling users to access their accounts effortlessly.




Doctor Dashboard – SwapDrill

A centralized interface where doctors can view and manage their assigned patients efficiently. Doctors can create new patient records, review existing ones, and reassign patients to other doctors with ease.



Assigning Patient to Another Doctor



Assigned Patient Details


```

Success: {success: true, message: 'Patient registered successfully', data: {}}
  data:
    clientId: "88FvGKSG5vqY90hXI2t9PUsKx6peJqHJ"
    clientSecret: "LuUeBtTfQNW--Bs_pM80p8rpz9m8zvznzLkqguISiSIXyN5y3Rwo3FRjVP0Iihix"
    email: "renu@gmail.com"
    firstName: "Renu"
    id: "cd80cda3-d3be-4161-8e00-e08d4461926b"
    lastName: "v"
    [[Prototype]]: Object
  message: "Patient registered successfully"
  success: true
  [[Prototype]]: Object


```

<D:\office\o-auth\mac-ntForm\page.tsx:121>

When creating Patient we get the clientId and ClientSecret using this we will send the request to the machine



Sending the request to the machine for the access .if it have the permission it will return the token



it have been return the token all this happens in the machine to machine auth


```

D:\office\o-auth>node sensordata.js
Original Sensor Data:
{
  {
    timestamp: '2025-05-15T15:52:07.894Z',
    temperature: 25.6,
    humidity: 82,
    light: 506
  },
  {
    timestamp: '2025-05-15T15:52:07.895Z',
    temperature: 26.7,
    humidity: 76,
    light: 473
  },
  {
    timestamp: '2025-05-15T15:52:07.895Z',
    temperature: 32.9,
    humidity: 84,
    light: 716
  }
}

Encrypted Sensor Data:
{
  {
    encrypted: 'b3908548816c2b917a6a4b08d8ae5a8b:ba547ee55295ffe23bfcba9ac7a7c2ffacdcc53aa7b7b5f61ed8f576554e3098402b890750a76ddb1f629cd511d9317b0fb18ab819cda5f2eb2f50a061ba1c3f623283d116b32650008de90b703fdf4b9d44b4ad2e7f40c56c0b96dc7b441084'
  },
  {
    encrypted: '34608ec0da296eca61682d0aef1586eb:8a346bedb735bc6b187581b52946fa9505ffd0e8792f70d5cd4941245bdf97945504103a4b010fefcd81a1557928e83964191e28b45c73f1ec2d70b5c1248cd5577b120e12cdedfeba69615f812454433a7e0a1df26bc75c1fb31f6fc7c8a0'
  },
  {
    encrypted: 'ec6bf7ba82d102c9c45f993e97ae0f6:06c2d116a524a47fe0510e1a47ebd20640286216deecd04bffe960f443d4643d1353ab66e33927c40a7464e27b664900ce8e9d53c40c4394ae4fc46bbdb2a0559baf30550b3a1ce6014eab7dedad6a9a1ae5dd1583d5e90daf0b3795449c64'
  }
}

```

Then we generate the Sensor data as the encrypted value



Send the value and decrypte and store it

HTTP <http://localhost:3000/api/sensore> Save No Environment

POST <http://localhost:3000/api/sensore> Send

Params Authorization Headers (11) Body Scripts Settings Code Cookies

☐ none ☐ form-data ☐ x-www-form-urlencoded ☒ raw ☐ binary ☐ GraphQL Text

```
1 {
2   "encryptedDataArray": [
3     ...
4   ]
5 }
```

Body Cookies Headers (6) Test Results Status: 200 OK Time: 8.67 s Size: 311 B

Pretty Raw Preview JSON

```
1 {
2   "message": "Sensor data saved successfully",
3   "savedRecords": 3
4 }
```

Patient x User x Close all

Filters None Fields All Showing 7 of 7 Add record

	id A	firstName A	lastName A	appointmentLocation A	email A	phoneNumber A	reasonForApp
☐	0d73ecb0-50c3-4f43-a0e6...	stanly	M	Fever	stanly@gmail.com	9363384134	Madurai
☐	4ca606e4-c47f-4cdf-9d2c...	rizwan	n	Fever	rizwan@gmail.com	9898983938	Coimbatore
☐	602c9f18-709d-4257-8b9d...	kumar	M	Fever	kumar@gmail.com	987654321	Madurai
☐	8ab34a36-75a0-4493-8016...	stephen	m	fever	stephen@gmail.com	7788996655	Trivandrum
☐	bb895524-9c2f-49da-b36a...	Rani	M	fever	rani@gmail.com	987654324	Madurai
☐	cd80cda3-d3be-4161-8e00...	Renu	C	Fever	renu@gmail.com	9363384128	3 SensorData
☐							
☐	e9cee210-03ed-46de-a8a1...	anzil	m	cold Fever	anzil@gmail.com	827383732732	Coimbatore

	id	timestamp	temperature	humidity	light	patientId
☐	search id...	search timestamp...	search temperature...	search humidity...	search light...	search patier
☒	c88a1eda-26db-4504-89c6...	2025-05-15T15:52:07.894Z	25.6	82	506	cd80cda3-d3
☒	93187bed-b2e0-4c9d-9195...	2025-05-15T15:52:07.895Z	26.7	76	473	cd80cda3-d3

Open in new tab Skip to unconnected records

This are the Decrypted data have been store to the patient